



MADE BY FARSEE



## Sealing Fit Nose Bridge



- Sealing-fit design can avoid air-leaking.
- No fog on glasses.
- No pollutant gets through.

## Mask Rope



- The ear lanyard has very good flexible and elastic characters.
- It can fit different sizes and wear comfortably.



Waterproof &  
Breathability



## NON-WOVEN FABRICS

Soft skin friendly, isolated particles

## HOT AIR COTTON


Antibacterial, odor and fog proof

## MELT BLOWN FABRIC

Electrostatic adsorption, bacteria filtration

## NON-WOVEN FABRICS

Skin friendly, breathable



**THICK NON-WOVEN FABRIC**

Filter over 95% of dust particles

**MELTBLOWN CLOTH**

Filter virus by efficient electrostatic adsorption

**PP NON-WOVEN ISOLATION LAYER**

Filter ommissive dust



**PP NON-WOVEN ISOLATION LAYER**

Filter micro-dust

**NON-WOVEN FABRIC**

Skin-friendly  
Isolate main particles






Great Breathability






MADE BY FARSEE



*Avoid Various Pollution*





***Suitable to ...***

*hospital, subway,  
office, elevator,  
Super market,  
school, etc.*




***OEM  
available***



# Wearing Indication



1, Wear the mask rope onto ears.



2, Press the nose bridge tightly.



3, Adjust the rope into a suitable position.




4, Pull the mask to fit cheek and chin.





## Head-band for choice



1. Open the mask.



2. Adjust the position.



3. Wear the mask by  
two head bands.



4. The lower band is  
on the neck.



5. The upper band is  
on the head.



6. Press the nose clip  
and adjust the mask to  
face.